



CardioGraphe™



CardioGraphe™



For first-line individualized decision-making in CAD

The only dedicated cardiovascular CT

The only compact system dedicated to cardiovascular CT that increases 1st line patients' access to tailored CAD management (prevention, diagnosis and treatment) because it is an easy-to-use, evidence-based tool to aid clinical decision-making. These elements provide easy patient and scan preparations to allow for an increase in challenging patient access by up to 15%.



Innovative design

CardioGraphe's™ novel technologies allow for a cardiac-centric patient experience. These include:

- **Stereo CT**—dual overlapping beams help reduce cone beam artifacts and enhance coverage at spatial resolutions of 0.28 mm
- **Focused FOV of 25cm**—provides 10~20 times dose reduction in areas outside of cardiac anatomy
- **Ultra-short geometry**—allows for faster rotation speeds of 0.24 seconds and more photon flux for a given tube power
- **Compact Design**—the CardioGraphe™ fits through standard doorways and elevators, and sits in rooms as small as 15m²

An enhanced cardiovascular care pathway

The CardioGraphe™ offers a multitude of advanced applications and workflow improvements to deliver a streamlined process for cardiac CT exams.

- **ASiR-CV**—reduces image noise and can be used to reduce dose while improving IQ with CardioGraphe's™ 25cm FOV and 140mm detector coverage
- **SmartScore 4.0**—advanced imaging software that detects, quantifies, and scores cardiac calcium plaque burden
- **CardIQ Xpress Reveal**—anatomical, functional, and perfusion information that aids in quick diagnosis
- **TAVI Analysis**—comprehensive work-up for TAVI procedure planning in a single workflow
- **CardEP**—streamlined post-processing for improved electrophysiology procedures





One Beat Cardiac Imaging

Motion-free cardiovascular CT at any heart rate


Challenge

- Patients with variable or high heart rhythms can necessitate the use of pharmaceuticals to stabilize cardiac motion
- Many of this patient base (about 15%) may have contraindications or arrhythmias which can produce suboptimal imaging even after following protocol measures


Solution

Cardiac structures move quickly and irregularly. Wide detector coverage of 140mm with rotation speeds up to 0.24 seconds help enable unlimited one-beat cardiac acquisitions


- Removes banding artifacts caused by multi-beat studies
- Provides an opportunity to expand your patient base



Unlimited one-beat cardiac, with wide coverage and fast rotation speeds, reduces the issues with multi-beat scanning by capturing the cardiac anatomy in one R-R cycle at any heart rate and patient size.




64/128 slice CT



- Can have limited IQ with stair and motion artifacts
- High and variable heart rates difficult to image

Unlimited 1-Beat GE Scanners



- Acquire and freeze motion in the heart in one rotation
- Simplify the cardiac workflow for technologists



SnapShot Freeze

Improves effective temporal resolution with intelligent motion correction

Challenge

- Cardiac motion can influence the accurate assessment of coronary arteries, especially in patients with high or variable heart rates, or degrade the image quality for cardiac evaluation.
- Technical developments in CT systems have improved the temporal resolution, including increased gantry rotation speeds, larger coverage, or dual-source CT approach. However residual cardiac motion remains. Additional motion correction algorithm is needed.

Solution

- Snapshot Freeze helps reduce motion blurring in vessels by up to a factor of six*.
- Characterizes motion within a single heart cycle, so it is not susceptible to beat-to-beat inconsistencies or heart period /gantry period resonance points which can plague multi-sector (i.e., multi-heart cycle) reconstruction.





Figure 2. Primary components of coronary motion correction.

Results

SnapShot Freeze provides automated cardiac vessel correction



* "Intelligent motion correction with SnapShot™ Freeze is designed to reduce blurring artifacts due to motion in coronary vessels that cannot be addressed by gantry speed alone. A 6x improvement of motion-blur reduction while maintaining high spatial resolution is demonstrated in cardiac phantom testing."



Building a world that works

© 2022 General Electric Company. All rights reserved. GE Healthcare is a leading global medical technology and digital solutions innovator. GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), and [Insights](#), or visit our website www.gehealthcare.com for more information.

JB19553XX