



GE HealthCare

Accessing your AW Server platform has never been better

AW Server's web client provides browser-based access to a refreshed and redesigned worklist interface and the ability to deploy the latest advanced visualization applications*



Access from anywhere

- ✔ Enables users to access the AW Server platform via commonly used web-browsers (requires internet connection and VPN login credentials)



Simplify maintenance with effortless updates

- ✔ No need to install thin-client application on users' desktops
- ✔ Updates or upgrades deployed automatically
- ✔ Consistent user experience



Enhanced interface







- ✔ Optimized worklist filtering
- ✔ Consistent user interface with other GE HealthCare products
- ✔ One-click review of relevant longitudinal patient studies

* Applications must be purchased separately from AW Server. Only available in certain countries based on regulatory approval



GE HealthCare

The minimum hardware configuration is as follows:

Processor 	Memory 	Disk drive 
2.33 GHz Pentium®4 minimum, Pentium dual core processors recommended for optimal performance	4 GB minimum	500 MB free space available
Network card 	Internet connection 	Mouse 
100 Mbps minimum (1000 Mbps recommended)	Customer-provided IPSEC VPN, for internet/WAN operation	Two or three-button mouse (three button suggested for best functionality)

Interested in learning more?

Contact your GE HealthCare sales representative.